

ABL800

BASIC

Meeting your essential needs



SIMPLY DEDICATED

RADIOMETER
COPENHAGEN



The right analyzer when measuring basic STAT parameters

The ABL800 BASIC is the ideal choice for you who are looking for efficiency and value in your critical care testing.

The ABL800 BASIC offers determination of the most important critical care parameters, including optimal values for total hemoglobin and oxygen saturation - with maximum accuracy and minimum effort.

Easy to use, the ABL800 BASIC provides you with fast and reliable results for fast patient treatment.

In addition, automatic quality control, automated sample aspiration and IT connectivity help reduce your workload and minimize errors.

pH pCO₂ pO₂ ctHb sO₂ cK⁺ cNa⁺ cCa²⁺ cCl⁻ cGlucose cLactate



Fewer steps and better workflow

Ease of use

- Automated sample aspiration
- Integrated barcode reader
- Inlet design reduces risk of air bubbles
- Configurable interface for custom profiles
- User-definable shortcuts
- Printouts based on user-defined values

Improve workflow

- User- access control
- Maintenance scheduler
- Stand-by option
- Parameter bar provides analyzer status at a glance

Save time and reduce errors

IT connectivity will further enhance the integration of blood gas in your processes and help reduce data management costs.

- Automatic storage of patient, QC, calibration and system data
- Two-way communication with HIS/LIS eliminates transcription errors
- Remote control and monitoring of analyzer status, e.g. through the RADIANCE analyzer data management system



Configurable interface for custom profiles.

IT connectivity

Computer specifications

Intel Celeron Processor
128 MB RAM
Hard disk 40 GB
CD/RW/DVD drive
TFT 10.4" VGA color touch screen
Dedicated 80386 CPU for wet-section operations

Data capacity

Patient results:	2000
Calibration results:	1000
QC results:	1500
System messages and service registrations:	5000

Communication

Access to Local Area Network for backup, etc.:
using PC network operating systems
supporting Windows®XP

Output protocols:

High-level protocols

ABL700-compatible ASTM (E1394-91)
ABL700-compatible HL7 (Version 2.2)
ABL5xx-compatible ASTM (E1394-91)

Low-level protocols


ASTM (E1381-95)
Radiometer network protocol (TCP/IP only)
Raw (serial only)

Transport layer

TCP/IP
RS232

RADIANCE

Communication via TCP/IP



Online assistance with audio-enhanced video tutorials, troubleshooting and help.

Maximize your sample success rate

The ABL800 BASIC offers an intuitive sampling process, where guides and narrated on-screen video tutorials offer immediate help for the user.

A wide on-screen choice of measuring modes ensures sample accuracy for volumes as low as 35 µL. Micromodes offer a fast response and high success rates for the precious capillary samples.

Mode	Sample volume (µL)	Measuring time (sec)	Cycle time (sec)
All parameters	195	80	150
All parameters, micro	95	135	200
pH + BG + Oxi, micro	85	80	145
pH + BG + Oxi, micro	55	100	170
Glu + Lac, micro	35	80	145
Oxi, micro	35	80	145

Mode	Sample volume (mL)	Measuring time (sec)	Cycle time (sec)
Expired air	15	65	170

Small sample volumes.

Save time with automatic quality control


On the ABL800 BASIC, quality control is completely automated, so you can schedule it according to your needs.

In addition, the ABL800 BASIC automatically saves QC data and displays Levey-Jenning's plots. All QC data is available electronically with easy search options.

For additional verification of analyzer performance, simply submit your QC data to Radiometer's online peer-group comparison service Worldwide DATACHECK.

AutoCheck automatic quality control

Number of ampoules in carousel:	0-20
Positioning of ampoules in carousel:	Random
Lot change:	2 lots of same level possible at the same time
Liquid volume in ampoule:	700 µL
Expiration of ampoules:	24 months at 25 °C/77 °F (including 15 days at up to 32 °C/90 °F)
Conditioning time (from room temperature):	15 min with filled carousel
Scanning time:	< 30 sec with filled carousel
Cycle time:	< 2 min 40 sec
Manual QC measurement possible:	Yes
Remote control:	Remote monitoring and start of measurement, e.g. via RADIANCE



Essential parameters

Measured parameters

On the ABL800 BASIC, you get all the essential parameters in one sample.

The high-performance oximeter on the ABL800 BASIC ensures optimal values for total hemoglobin (tHb) and oxygen saturation, thus providing you with a clear picture of the patient's levels of oxygen uptake, transport and release.

	Type	Parameters	Units	Measuring ranges
Default parameters	pH	pH cH ⁺	pH scale nmol/L	6.300-8.000 10.0-501
	Blood gas	pCO ₂	mmHg	5.0-250
			kPa	0.67-33.3
			Torr	5.0-250
Additional parameters	Electrolyte*	cCl ⁻	mmol/L	7-350
			meq/L	7-350
		cCa ²⁺	mmol/L	0.20-9.99
			meq/L	0.40-19.98
	Metabolite*	cK ⁺	mmol/L	0.5-25.0
			meq/L	0.5-25.0
		cNa ⁺	mmol/L	7-350
			meq/L	7-350
	Metabolite*	cGlu	mmol/L	0.0-60
			mg/dL	0-1081
		cLac	mmol/L	0.0-30
			mg/dL	0-270

* The total number of electrolytes and metabolites is limited to four.



Derived parameters

Type	Definition
pH(T)	pH of blood at patient temperature
pCO ₂ (T)	Carbon dioxide tension of blood at patient temperature
cHCO ₃ ⁻ (P)	Concentration of hydrogen carbonate in plasma
cBase(B)	Concentration of titrable base of blood (actual base excess)
cBase(B,ox)	Actual base excess at 100 % oxygen saturation
cBase(Ecf)	Concentration of titrable base of extracellular fluid (standard base excess)
cBase(Ecf,ox)	Standard base excess at 100 % oxygen saturation
cHCO ₃ ⁻ (P,st)	Concentration of hydrogen carbonate in plasma of standardized blood (standard bicarbonate)
cH ⁺	Concentration of hydrogen ions in blood
cH ⁺ (T)	Concentration of hydrogen ions in blood at patient temperature
ctCO ₂ (P)	Concentration of total carbon dioxide in plasma
ctCO ₂ (B)	Concentration of total carbon dioxide of whole blood (CO ₂ content)
pH(st)	pH of standardized blood (pCO ₂ = 40 mmHg)
pO ₂ (T)	Oxygen tension of blood at patient temperature
pO ₂ (A)	Oxygen tension of alveolar air
pO ₂ (A,T)	Oxygen tension of alveolar air at patient temperature
p50	Oxygen tension at 50 % saturation of blood
p50(T)	Oxygen tension at 50 % saturation of blood at patient temperature
p50(st)	Oxygen tension at 50 % saturation of blood at standard conditions for pH, pCO ₂ , FCOHb, FMetHb, FHbF at 37 °C
pO ₂ (A-a)	Difference of oxygen tension of alveolar air and arterial blood
pO ₂ (A-a,T)	Difference of oxygen tension of alveolar air and arterial blood at patient temperature
pO ₂ (a/A)	Ratio of oxygen tension of arterial blood and alveolar air
pO ₂ (a/A,T)	Ratio of oxygen tension of arterial blood and alveolar air at patient temperature
pO ₂ (a)/FO ₂ (I)	Oxygen tension ratio of arterial blood to the fraction of oxygen in inspired air
pO ₂ (a,T)/FO ₂ (I)	Oxygen tension ratio of arterial blood at patient temperature to the fraction of oxygen in inspired air
cCa ²⁺ (pH=7.40)	Concentration of ionized calcium in plasma at pH 7.40
Anion Gap(K ⁺)	Concentration difference of K ⁺ + Na ⁺ and Cl ⁻ + HCO ₃ ⁻
Anion Gap	Concentration difference of Na ⁺ and Cl ⁻ + HCO ₃ ⁻
DO ₂	Oxygen delivery
Hct	Fraction of the volume of erythrocytes in the volume of whole blood
pO ₂ (x)	Oxygen extraction tension of arterial blood
pO ₂ (x,T)	Oxygen extraction tension of arterial blood at patient temperature
ctO ₂ (B)	Total oxygen concentration of blood (O ₂ content)
ctO ₂ (a-v)	Total oxygen concentration difference between arterial and mixed venous blood
BO ₂	Oxygen capacity of hemoglobin. The maximum concentration of oxygen bound to hemoglobin in blood, saturated so that all deoxyhemoglobin is converted to oxyhemoglobin
ctO ₂ (x)	Extractable oxygen concentration of arterial blood
FShunt	Volume fraction of shunted venous blood in arterial blood
FShunt(T)	FShunt at patient temperature
RI	Respiratory Index
RI(T)	Respiratory Index at patient temperature
VO ₂	Oxygen consumption
mOsm	Plasma osmolality
Qx	Oxygen compensation factor of arterial blood
Q _t	Cardiac output
V(B)	Volume of blood
sO ₂	Saturation
FO ₂ Hb	Fraction of oxyhemoglobin in total hemoglobin in blood

Extendability with unlimited input parameters

Type	Definition	Type	Definition
Patient ID	Patient identification number	Patient name	Name of the patient
Patient height	The height of the patient	Physician	Name of the physician
Patient department	Which department the patient is from	Operator	Name of the operator
T	Patient temperature	Operator department	Which department the operator is from
Sample type	Arterial, venous, mixed venous, capillary, prof. test, cal. verification	p50(st)	Oxygen tension at 50 % saturation of blood at standard conditions for pH, pCO_2 , FCOHb, FMetHb, FHbF at 37 °C
Patient note	Notes about the patient or sample	RQ	Respiratory quotient
Patient weight	The weight of the patient	FO ₂ (I)	Fraction of oxygen in dry inspired air
Patient Accession No.	Specific sample number	Q _t	Cardiac output
Patient age	Date of birth	VO ₂	Oxygen consumption
Patient sex	Male or female	VCO	Volume of carbon monoxide, input value for measurement of V(B)
Draw time	When the sample was taken	sO ₂ (v̄)	Oxygen saturation of hemoglobin in mixed venous blood
Date of birth	Patient date of birth	pO ₂ (v̄)	Oxygen tension of mixed venous blood
Sample site	Not specified, brachial left/right, femoral left/right, radial left/right, finger left/right, heel left/right, umbilical cord	ctHb	Total hemoglobin concentration (if not measured)
Patient birth weight	The weight of the newborn	FCOHb(1)	Used for determining blood volume
Patient gestational age	Period of intrauterine fetal development from conception to birth	FCOHb(2)	Used for determining blood volume
		Parameters can be set for mandatory input and are user-definable.	

Maximum uptime

Calibration data			
Automatic:	Default interval:	Interval options:	
1-point cal.	4 hours	after measurement, 30 min, 1, 2, 4 hours	
2-point cal.	4 hours	after measurement, 1, 2, 4 hours	
System alignment	24 hours		
Cleaning	24 hours	8, 24 hours	
Manual:			
tHb calibration	3 months	never, 7 days, 1, 2, 3, 4, 6 months	



Adjustable color touch screen.




Integrated barcode reader.



Automatic sample aspiration.

For more information about Radiometer, our products and services, visit www.radiometer.com.

Read how hospitals around the world are optimizing their blood gas testing. Visit Radiometer's knowledge site www.bloodgas.org.



Additional information



EMC Emission

Indicates compliance with the IVD Directive 98/79/EC

EMC Immunity

Patents

Approvals
Languages

The equipment complies with the emission requirements for Class B equipment in EN 61326-1:

Electrical equipment for measurement, control and laboratory use
- EMC requirements part 1: General requirements

The equipment complies with the immunity requirements in EN 61326-1:

Electrical equipment for measurement, control and laboratory use
- EMC requirements part 1: General requirements

One or several of the following patents and patent applications may apply:

European Patent Nos.: EP210417, EP212126, EP889951, EP944731

Japanese Patent Nos.: JP2070313, JP2972351, JP3285879, JP3369547

German Patent Nos.: DE3673910, DE3686855, DE69729185

Austrian Patent No.: AT56271

Danish Patent No.: DK155764

European Patent Application Nos.: EP1084398, EP1086366, EP1273920

Japanese Patent Application Nos.: JP2000/507457, JP2002/518670

International Patent Application Nos.: WO2004/057304, WO2005/052596

Other patents pending.

In compliance with IEC 61010-1 Installation Category II

English, German, Spanish, Italian, French, Chinese, Portuguese, Danish, Greek, Hungarian, Estonian, Russian and Lithuanian.

Dimensions

Width	70 cm	28 in
Height	55 cm	22 in
Depth	55 cm	21 in
Weight	33.9 kg	75.5 lbs
Warm-up time	Cold start: 25 min typical. Warm start: 5 min	
Ambient temperature	15-32 °C/59-90 °F	
Relative humidity	20-80 %	
Thermostatting	pH and blood gases, 37.0 °C ± 0.15 °C/98.6 °F ± 0.3 °F	
	Electrolytes and metabolites, 37.0 °C ± 0.25 °C/98.6 °F ± 0.5 °F	
Spectrometer	128-wavelength measurement	
Hemolyzer frequency	30 KHz intracuvette hemolysis	
Barometer	450-800 mmHg	
Power	100-240 VAC, 50-60 Hz, 250 VA	

IVD

Sales companies:

Country:	Radiometer representative:	Telephone:
Australia:	Radiometer Pacific Pty. Ltd.	(+61) 3 9259 2222
China:	Radiometer China	(+86) 21 6354 3387
Denmark:	Radiometer Danmark	(+45) 38 27 28 29
France:	Radiometer S.A.S.	(+33) 1 49 44 35 50
Germany:	Radiometer GmbH	(+49) 2154 818-0
Ireland:	Radiometer Ireland Ltd.	(+353) 1 888 3611
The Netherlands:	Radiometer Nederland BV	(+31) 79 361 4593
New Zealand:	Radiometer New Zealand	(+64) 9 574 1400
Poland:	Radiometer Sp.z.o.o.	(+48) 22 518 02 40
Portugal:	Radiometer Ibérica, S.A.	(+351) 214 12 39 70
Spain:	Radiometer Ibérica, S.A.	(+34) 91 655 99 50
Switzerland:	Radiometer GmbH	(+41) 044 723 38 60
United Kingdom:	Radiometer Ltd.	(+44) 1293 517 599
Other countries:	Radiometer International Sales Division	(+45) 38 27 38 27

E-mail:
sales@radiometer.com.au
pol@radiometer.com.cn
rdan@radiometer.dk
marketing@radiometer.fr
info@radiometer.de
sales@radiometer.ie
info@radiometer.nl
sales@radiometer.co.nz
info@radiometer.pl
rapor@radiometer.es
resp@radiometer.es
info@radiometer.ch
sales@radiometer.co.uk
isd@radiometer.dk

Website:
www.radiometer.com.au
www.radiometer.cn
www.radiometer.dk
www.radiometer.fr
www.radiometer.de
www.radiometer.ie
www.radiometer.nl
www.radiometer.co.nz
www.radiometer.pl
www.radiometer.pt
www.radiometer.es
www.radiometer.ch
www.radiometer.co.uk
www.radiometer.com

ABL800 BASIC is not available for sale in the USA or Canada

Data subject to change without notice.

Radiometer™, Red System™, ABL™, AutoCheck™ and RADIANCE™ are trademarks of Radiometer Medical ApS, Denmark. Windows® is a registered trademark of Microsoft Corporation.

